



REF

ES-09807

8

Fr.

10

cm sheath length

.035

inch dia. spring-wire guide

PSI Kit for use with 7 - 7.5 Fr. Catheters

Contents:	
1: Arrow-Flex® Sheath: 8 Fr. (2.5 mm tip ID) x 10 cm and Tissue Dilator	2: Towel
1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer	1: Scalpel: #11
1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle	5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
1: Injection Needle: 25 Ga. x 1" (2.54 cm)	1: Suture: 1-0 Silk with Straight Needle
1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW	1A registered trademark of CareFusion or one of its subsidiaries.
1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)	
1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe	
1: Syringe: 3 mL Luer-Lock	
1: Syringe: 5 mL Luer-Slip	
1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep	
1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ¹ One-Step Solution with Hi-Lite Orange™ Tint	
1: Obturator: 8 Fr.	
1: Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape	
1: SharpsAway® Disposal Cup	
1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration, with adhesive	

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
California Prop. 65
⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov
Caution: This Product Contains Dry Natural Rubber.
Store between 20 – 25°C (68 – 77°F).
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.



!

MD

i

⚠

⌘

2

STERILIZE

STERILE EO

⬭

☀

☔

⊠

LATEX

20°C (68°F) 25°C (77°F)

LOT

⏳

🏭

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058320 R01 (2022-02)



REF

ES-09807

8

Fr.

10

cm sheath length

.035

inch dia. spring-wire guide

PSI Kit for use with 7 - 7.5 Fr. Catheters