



REF NS-04100-S

20^{Ga.} 12^{cm} catheter length .025^{inch dia.} spring-wire guide

Arterial Catheterization

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP
 - 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
 - 1: Catheter: 20 Ga. x 1-1/4" (3.2 cm) Jelco®¹ Catheter over 22 Ga. Introducer Needle
 - 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
 - 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
 - 2: Syringe: 3 mL Luer-Lock
 - 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
 - 1: Prep Sponge Swab
 - 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- ¹A registered trademark of Medex, Inc.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex
www.teleflex.com

LBL059321 R01 (2022-02)



REF NS-04100-S

20^{Ga.} 12^{cm} catheter length .025^{inch dia.} spring-wire guide

Arterial Catheterization