

20 Ga.

3.81 cm
catheter
length

.025 inch dia.
spring-wire
guide

Contents:

- | | |
|--|---|
| 1: Indwelling Catheter: 20 Ga. x 3.81 cm Radiopaque Polyurethane with Integral Suture Wing | 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends) | 3: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Injection Needle: SafetyGlide™ 18 Ga. x 1-1/2" (3.81 cm) | 1: Dressing: Tegaderm® ³ I.V. 11.5 cm x 8.9 cm |
| 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW | 1: Surgical Apparel: Impervious Gown |
| 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm) | 1: Surgical Apparel: Bouffant Cap |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Syringe: 3 mL Luer-Lock | 1: Suture: 3-0 Silk with Curved Needle |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | ¹ A trademark of Becton, Dickinson and Company. |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ² One-Step Solution | ² A registered trademark of CareFusion or one of its subsidiaries. |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ² One-Step Solution with Hi-Lite Orange™ Tint | ³ A registered trademark of 3M Company. |
| 1: SharpsAway® Disposal Cup | |
| 1: SharpsAway® II Locking Disposal Cup | |
| 1: Drape: 41" x 55" with 4" fenestration | |
| 1: Needle Holder | |
| 1: Forceps | |
| 1: Safety Scalpel #11 | |

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).


Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL059358 R01 (2023-05)

REF ASK-04120-HF

20 Ga.

3.81 cm
catheter
length

.025 inch dia.
spring-wire
guide

ARROW®