



REF PW-18045

**.018** inch dia.  
spring-wire  
guide


**45** cm  
spring-wire  
guide

**Spring-Wire Guide .018" (0.46 mm) diameter x 17-3/4" (45 cm) wire / Proximal End: Solid Wire .018" diameter / Distal End: Coil .018" diameter x 9 cm**

Contents:


1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

Rx only [www.teleflex.com/IFU](http://www.teleflex.com/IFU)



 EU Authorized Representative and Importer:

 Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

 Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL060572 R01 (2023-04)

LOT



REF PW-18045

**.018** inch dia.  
spring-wire  
guide

**45** cm  
spring-wire  
guide