



REF HO-14702

7 Fr. 2 Lumen

30 cm catheter length

23 cm insertable catheter length with Cath-Gard®

HANDS-OFF® Two-Lumen CVC for use only with Arrow One-Piece PSI System

Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm, Cath-Gard® Catheter Contamination Shield
- 2: Dust Cap: Non-Vented

All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



CAS # 7440-48-4



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.57	1917
Proximal (16 Ga.)	0.62	2000

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL061678 R01 (2023-08)

LOT



7 Fr. 2 Lumen

30 cm catheter length

23 cm insertable catheter length with Cath-Gard®

REF HO-14702

HANDS-OFF® Two-Lumen CVC for use only with Arrow One-Piece PSI System