

REF EU-00085-CC

8.5 Fr. Catheter Clamp with Fastener

Contents:

- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU


Not made with natural rubber latex. Fluid path components are non-pyrogenic.

en fr de it pt es sv



 EU Authorized Representative and Importer:

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
2797

 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

LOT



LBL062808 R00 (2021-12)

Teleflex®
www.teleflex.com