



REF CK-04018

ARROW®

18^{Ga.} 12^{cm} catheter length .025^{inch dia.} spring-wire guide

Arterial Catheterization Set

Contents:

- 1: Catheter: 18 Ga. x 12 cm Radiopaque Polyurethane with Integral Suture Wing
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.



LOT



EU Authorized
Representative and
Importer:

EC REP

CE 2797

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Arrow International LLC

Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL062848 R00 (2022-01)