



REF GH-04150

20^{Ga.} 20^{cm} catheter length .025^{inch dia.} spring-wire guide

Femoral Artery Catheterization Set

Contents:

- 1: Catheter: 20 Ga. x 20 cm Radiopaque Polyurethane
 - 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
 - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/ IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.



EU Authorized Representative and Importer:

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LOT



L.B.L062866 R00 (2022-01)



20^{Ga.} 20^{cm} catheter length .025^{inch dia.} spring-wire guide

REF GH-04150

Femoral Artery Catheterization Set