



REF RC-09700

7^{Fr.} 5.08^{cm} catheter length .025^{inch dia.} spring-wire guide

RIC[®] Rapid Infusion Catheter Exchange Set

Contents:

- 1: Sheath: 7 Fr. x 2" (5.08 cm) Arrow-Flex[®] Radiopaque Polyurethane with Tissue Dilator
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: Scalpel: #11

All components are CE 2797 unless otherwise noted.



Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



EU Authorized Representative and Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex[®]
www.teleflex.com

LBL064093 R00 (2022-02)

LOT

