



REF ES-04522-BZ

22 Ga. 1 Lumen 10 cm catheter length .021 inch dia. spring-wire guide

# One-Lumen CVC

**Contents:**

- 1: One-Lumen Catheter: 22 Ga. (0.9 mm OD) x 10 cm
- 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 4 Fr. (1.35 mm) x 7.6 cm
- 1: Tissue Dilator: 4 Fr. (1.35 mm) x 4.1 cm
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**

**Not made with natural rubber latex.**  
Fluid path components are non-pyrogenic.



en pt

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (22 Ga.)	0.20	568

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL064467 R00 (2022-03)

LOT



22 Ga. 1 Lumen 10 cm catheter length .021 inch dia. spring-wire guide

REF ES-04522-BZ

# One-Lumen CVC