

Product of xxx  
Packaged in xxx



REF RHY-177-TTSAP (IPNXXXXXX)



# VPS TipTracker™ Stylet Accessory

For exclusive use with VPS Rhythm™ Device. Consult Operator's Manual for Use.

### Contents:

- 1: VPS TipTracker™ Stylet with T-Port Connector: 0.015" (0.38 mm) dia. x 29" (74 cm)
- 1: Foam Electrodes (3 per pack)
- 1: Equipment Cover: 13 cm x 244 cm (5" x 96")
- 1: T-Piece Cover

All components are CE 2797 unless otherwise noted.

**Rx only Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

California Prop. 65

**WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

EU Authorized Representative and Importer:

Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
2797 Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**teleflex**  
[www.teleflex.com](http://www.teleflex.com)

LBL067352 R01 (2024-05)