



REF EDC-00622-M (IPNXXXXXX)

22 Ga. catheter 6 cm catheter length .010 inch dia. wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

Contents:

- 1: Catheter+ System: 22 Ga. (0.91 mm) x 60 mm (2.36") Arrow™ Endurance™ over 24 Ga. (0.60 mm) TW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
  - 1: 10 mL Luer-Lock Pre-Filled Saline Syringe
  - 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint
  - 1: Label: Venous/Arterial Placement
  - 1: Tourniquet
  - 2: Gauze Pad: 4" x 4" (10 cm x 10 cm)
  - 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm
  - 1: Transducer Cover with Gel and Rubber Bands
- <sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
California Prop. 65  
⚠️ WARNING: Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
\*\*\*MR Safe  
Not made with natural rubber latex.  
Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).  
Fluid path components are non-pyrogenic.





LOT Lot Number

YYYY-MM-DD

UDI

YYYY-MM-DD

(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

CAS # 7440-48-4 MR



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)  
LBL069108 R02 (2025-02)



REF EDC-00622-M

22 Ga. catheter 6 cm catheter length .010 inch dia. wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

