



REF EDC-00822-M

(IPNXXXXXX)

22 Ga. catheter

8 cm catheter length

.010 inch dia.
wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

Contents:

- 1: Catheter,††† System: 22 Ga. (0.91 mm) x 80 mm (3.15") Arrow™ Endurance™ over 24 Ga. (0.60 mm) TW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- 1: 10 mL Luer-Lock Pre-Filled Saline Syringe
- 1: 3 mL Applicator 2% CHG and 70% IPA ChlorPrep®† One-Step Solution with Hi-Lite Orange™ Tint
- 1: Label: Venous/Arterial Placement
- 1: Tourniquet
- 2: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm
- 1: Transducer Cover with Gel and Rubber Bands

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65


WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

###MR Safe

Not made with natural rubber latex.

Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.



 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

LOT Lot Number

YYYY-MM-DD

UDI



(01)XXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



www.teleflex.com/IFU

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LBL069111 R02 (2025-02)



REF EDC-00822-M

22 Ga. catheter

8 cm catheter length

.010 inch dia.
wire guide



Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

