



REF EDC-00822-HC

(IPNXXXXXX)

22 Ga. catheter

8 cm catheter length

.010 inch dia. wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

Contents:

- 1: Catheter System: 22 Ga. (0.91 mm) x 80 mm (3.15") Arrow™ Endurance™ over 24 Ga. (0.60 mm) TW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
 - 1: 10 mL Luer-Lock Pre-Filled Saline Syringe
 - 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint
 - 1: Arm Drape: 36" x 46" with 3" x 5" Fenestration
 - 1: Label: Venous/Arterial Placement
 - 1: Paper Tape Measure
 - 1: Tourniquet
 - 2: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Dressing: Tegaderm® CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
 - 1: Transducer Cover with Gel and Rubber Bands
 - 1: Dressing: STATLOCK® PICC Plus Catheter Stabilization Device
 - 1: Surgical Apparel: Bouffant Cap
 - 1: Surgical Apparel: Mask
- ¹A registered trademark of CareFusion or one of its subsidiaries.
²A registered trademark of 3M Company.
³A registered trademark of C. R. Bard, Inc.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
California Prop. 65
⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov
***MR Safe
Not made with natural rubber latex.
Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).
Fluid path components are non-pyrogenic.



Lot Number

Lot Number

UDI

(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

CAS # 7440-48-4



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com
LBL069116 R02 (2025-02)



www.teleflex.com/IFU



REF EDC-00822-HC

22 Ga. catheter

8 cm catheter length

.010 inch dia. wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

