



REF EDC-00620-HB

(IPNXXXXXX)

20 Ga. catheter

6 cm catheter length

.010 inch dia. wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

Contents:

- 1: Catheter*** System: 20 Ga. (1.10 mm) x 60 mm (2.36") Arrow™ Endurance™ over 22 Ga. (0.8 mm) RW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- 1: 10 mL Luer-Lock Pre-Filled Saline Syringe
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution with Hi-Lite Orange™ Tint
- 1: Arm Drape: 36" x 46" with 3" x 5" Fenestration
- 1: Label: Venous/Arterial Placement
- 1: Paper Tape Measure
- 1: Tourniquet
- 2: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm
- 1: Dressing: BIOPATCH®²
- 1: Transducer Cover with Gel and Rubber Bands
- 1: Dressing: STATLOCK®³ PICC Plus Catheter Stabilization Device
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask

¹A registered trademark of CareFusion or one of its subsidiaries.
²A registered trademark of Johnson & Johnson Corporation.
³A registered trademark of C. R. Bard, Inc.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
California Prop. 65
⚠️ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov
***MR Safe
Not made with natural rubber latex.
Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).
Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

UDI

YYYY-MM-DD

(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

CAS # 7440-48-4

MR



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com
LBL069117 R02 (2025-02)



www.teleflex.com/IFU



REF EDC-00620-HB

20 Ga. catheter

6 cm catheter length

.010 inch dia. wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

