



REF EDC-00818-HB

(IPNXXXXXX)

18 Ga. catheter

8 cm catheter length

.010 inch dia. wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

Contents:

- 1: Catheter System: 18 Ga. (1.30 mm) x 80 mm (3.15") Arrow™ Endurance™ over 20 Ga. (1.0 mm) HW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- 1: 10 mL Luer-Lock Pre-Filled Saline Syringe
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 1: Arm Drape: 36" x 46" with 3" x 5" Fenestration
- 1: Label: Venous/Arterial Placement
- 1: Paper Tape Measure
- 1: Tourniquet
- 2: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm
- 1: Dressing: BIOPATCH®<sup>2</sup>
- 1: Transducer Cover with Gel and Rubber Bands
- 1: Dressing: STATLOCK®<sup>3</sup> PICC Plus Catheter Stabilization Device
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.  
<sup>2</sup>A registered trademark of Johnson & Johnson Corporation.  
<sup>3</sup>A registered trademark of C. R. Bard, Inc.

Rx only  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
California Prop. 65  
⚠️ WARNING: Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
\*\*\*MR Safe  
Not made with natural rubber latex.  
Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).  
Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

UDI

YYYY-MM-DD

(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

CAS # 7440-48-4 MR



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)  
LBL069119 R02 (2025-02)



[www.teleflex.com/IFU](http://www.teleflex.com/IFU)



REF EDC-00818-HB

18 Ga. catheter

8 cm catheter length

.010 inch dia. wire guide

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System

