

ARROW™ EZ-IO™

INTRAOSSSEOUS VASCULAR ACCESS



Intraosseous Access Tray + Single-Use Power Driver



45 mm 15 ga

**45 mm x 1.8 mm
40 kg or over**

REF 9079T-VC-005
(IPNXXXXXX)

Contents:

- 1: EZ-IO™ Power Driver (Single-Use)
 - 1: EZ-IO™ 45mm 15 Ga Intraosseous Needle‡ Set
 - 1: 10 mL Luer-Lock Pre-Filled Saline Syringe
 - 1: 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Frepp® Solution
 - 1: SharpsAway® II Locking Disposal Cup
 - 1: EZ-IO™ Wristband
 - 1: EZ-Stabilizer™ Dressing
 - 1: EZ-Connect™ Extension Set (Priming Volume - Approx. 1.0 ml)
- ¹A registered trademark of CareFusion or one of its subsidiaries.

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex™
ArrowEZIO.com
24-Hour Clinical Support: 1.888.413.3104
Customer Service: 1.866.479.8500

Product of xxx
Packaged in xxx

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65
⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

‡MR Conditional
Not made with natural rubber latex.
Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).

Rx only

LOT Lot Number
 YYYY-MM-DD
 YYYY-MM-DD

(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL069292 R02 (2024-01)

Peel Here

ARROW™ EZ-IO™
INTRAOSSSEOUS VASCULAR ACCESS

REF 9079T-VC-005

Intraosseous Access Tray
+ Single-Use Power Driver

45 mm 15 ga
45 mm x 1.8 mm
40 kg or over