



REF CAN-04020-NSH (IPNXXXXXX)

20 Ga.

4.45 cm catheter length

.018 inch dia. spring-wire guide

Arterial Catheterization Kit

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
 - 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm)
 - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
 - 1: Syringe: 3 mL Luer-Lock
 - 1: SharpsAway® Disposal Cup
 - 1: SharpsAway® II Locking Disposal Cup
 - 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
 - 1: Towel: 13" x 18" (33 cm x 46 cm)
 - 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 - 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Dressing: I.V. Advanced 10 cm x 12 cm
 - 1: Transducer Cover: 14 cm x 147 cm
 - 1: Dressing: STATLOCK®² Arterial Select Catheter Stabilization Device 8.9 cm x 3.8 cm
 - 1: Surgical Apparel: Bouffant Cap
 - 1: Surgical Apparel: Mask with Eye Shield
- ¹A trademark of Becton, Dickinson and Company.
²A registered trademark of C. R. Bard, Inc.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of XXX
Packaged in XXX

Teleflex®
www.teleflex.com

LBL071764 R00 (2024-04)



REF CAN-04020-NSH

20 Ga.

4.45 cm catheter length

.018 inch dia. spring-wire guide



Arterial Catheterization Kit

