

Seal/Weld Area

REF RA-04018 (IPNXXXXXXXX)

Radial Artery Catheterization Set

Contents:

- 1: Catheter: 18 Ga. x 4.45 cm Radiopaque FEP over 20 Ga. TW Introducer Needle with Integral .025" (0.64 mm) dia. Spring-Wire Guide
- 1: Wing Clip: Suture

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

ARROW

STERILE EO



LATEX



CAS # 7440-48-4



MD



Teleflex[®]
www.teleflex.com

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

LBL074025 R00 (2024-01)

Product of xxx
Packaged in xxx

LOT LotNumber

Hourglass icon YYYY-MM-DD

Wavy line icon YYYY-MM-DD

(01)XXXXXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber



T-00053-100B Dimensions: 11.890" x 1.7549" Step: 7 across web @ 1.7549" Step 2 direction of web @ 11.890"

T-00053-100B

Dimensions: Width: 11.890" Height: 1.558"