

Seal/Weld Area

REF GH-04125 (IPNXXXXXX)

Arterial Catheterization Set

Contents:

- 1: Catheter: 20 Ga. x 3-1/16" (7.78 cm) with Integral Suture Wing
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip Both Ends)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

ARROW

STERILE EO



LATEX



MD



STERILE



Teleflex[®]
www.teleflex.com

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

LBL074030 R00 (2024-01)

Product of xxx
Packaged in xxx

LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD

(01)XXXXXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber



T-00053-100B Dimensions: 11.890" x 1.7549" Step: 7 across web @ 1.7549" Step 2 direction of web @ 11.890"

T-00053-100B

Dimensions: Width: 11.890" Height: 1.558"