



REF ASK-09803-FM1



8.5^{Fr.}

10^{cm}
sheath
length

.035^{inch dia.}
spring-wire
guide

PSI Kit for use with 7.5 - 8 Fr. Catheters

Contents:

- 1: Single-Lumen Infusion Catheter, SLIC®: 7 Fr. x 16 cm
- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Injection Needle: Eclipse™¹ 23 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm)
- 1: Syringe: 3 mL Luer-Lock
- 2: Syringe: 5 mL Luer-Slip
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: SharpsAway® Disposal Cup
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Safety Scalpel: #11
- 1: Dressing: Tegaderm®² 4" x 4-3/4" (10 cm x 12 cm)
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Transducer Cover: 14 cm x 147 cm
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of 3M Company.

Not made with natural rubber latex.

