



REF

CV-15123-F

(IPNXXXXXX)

12^{Fr.}

3^{Lumen}

20^{cm catheter length}

.035^{inch dia. spring-wire guide}

Three-Lumen CVC

Contents:

1:

Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 20 cm

Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXX

1:

Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXX

1:

Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW

Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXX

1:

Syringe: 5 mL Luer-Slip

Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

1:

Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm

Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXX

3:

Dust Cap: Non-Vented

Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](#)
Not made with natural rubber latex.
Fluid path components are non-pyrogenic.

- en

be

bg

hr

cs

da

nl

et

fi

fr

de

el

hu

is

it

lv

lt

no

pl

pt

ro

ru

sr

sk

sl

es

sv

tr

uk

CAS # 7440-48-4

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.4	2224
Medial (12 Ga.)	1.6	10616
Proximal (12 Ga.)	1.5	12168

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EC

REP

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx

www.teleflex.com

LBL076175 R00 (2024-10)

LOT

Lot Number

XXXX-MM-DD

XXXX-MM-DD

(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number