

Dimensions:
Width: 7 3/8"
Length: 12"

PB-04450-202

T-04450-202 Dimensions: 7.360" x 11.219"



REF PW-18045

(IPNXXXXXX)

.018

inch dia.
spring-wire
guide

45

cm
spring-wire
guide

Spring-Wire Guide .018" (0.46 mm) diameter x 17-3/4" (45 cm) wire / Proximal End: Solid Wire .018" diameter / Distal End: Coil .018" diameter x 9 cm

Contents:

1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

Rx only www.teleflex.com/IFU



CAS # 7440-48-4

en

fr

de

it

pt

es

sv



EU Authorized
Representative and
Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LBL076666 R00 (2026-09)

LOT

Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



REF PW-18045

.018

inch dia.
spring-wire
guide

45

cm
spring-wire
guide



Spring-Wire Guide .018" (0.46 mm) diameter x 17-3/4" (45 cm) wire / Proximal End: Solid Wire .018" diameter / Distal End: Coil .018" diameter x 9 cm