



REF SI-09870-E (IPNXXXXXX)



8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

### Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
  - 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
  - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
  - 1: Pressure Transduction Probe
  - 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
  - 1: Obturator: 8 Fr.
  - 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
  - 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- All components are CE 2797 unless otherwise noted.

**Rx only**  
**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
**Caution:** This Product Contains Dry Natural Rubber.  
 Fluid path components are non-pyrogenic.

en bg cs da nl fi fr de el hu it no pl pt  
 ro ru sk sl es sv tr

EU Authorized Representative and Importer:  
 Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Product of xxx  
 Packaged in xxx

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)  
 LBL076690 R01 (2025-07)



LOT Lot Number  
 YYYY-MM-DD YYYY-MM-DD

(01)XXXXXXXXXXXXXX  
 (17)YYMMDD  
 (11)YYMMDD  
 (10)Lot Number



REF SI-09870-E

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide



# Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters